


Container Card



‘Label’ your sterile container...

Container cards are an important aid to label and mark the container. Not only the content but also details like sterilization date, Expiry date, the person who packed the container can be written down on the card. The chemical process indicator is giving the indication whether it has been sterilized or not. The container card can be used for keeping records at the operating room when it is stored in the patient record.

REF: 903.101.0250

Container Card

Unit: 250 pieces

Characteristics:

- ◆ Dimensions: 68 x 37,5 mm.
- ◆ Suitable for steam sterilization processes.
- ◆ Can be fully recycled.
- ◆ Steam indicator conform to ISO 11140 - Type 1
- ◆ Clear color transition of the indicator

Other types and dimensions available upon request.

