

3-Line Label guns & Labels



Three lines ... the all information label

The Sterintech™ three line label gun is reliable, easy to use and durable. It will serve for years of intensive usage. Sterintech™ double labels for three line label guns are available with all type of sterilization indicators as well as a variety of colors. Colored labels are used to assist FIFO systems as well as easy recognition of sterile goods for a certain department.

Often the combination of a double label gun label and normal tape is a more cost-effective solution than steam tapes only. Next to that they offer better readable data and prevent writing errors.

Characteristics:

- ◆ Process Indicators (Steam, FO, EtO, Plasma) conform to ISO 11140 - 1 Type 1
- ◆ Clear transition of indicator into its end-color.
- ◆ Easy loading of labels into gun
- ◆ Double labels serve for an easy traceability and documentation
- ◆ Non-Toxic / Lead Free indicators
- ◆ Latex free label materials
- ◆ Works with identical labels as the Blitz label gun

For reference numbers and types please refer to the next side.



3-Line Label guns & Labels



Three lines ... the all information label

Sterintech™ Three line label gun:

REF: 603.203.0001	Three line label gun	Unit: 1 piece
REF: 609.104.0001	Ink roller for Three line label gun	Unit: 1 piece

Sterintech™ double label gun labels are standard available as follows. Bespoke products also for other label gun systems are available upon request:

REF: 606.010.7500	Label with Steam indicator	Unit: 7.500 labels
REF: 606.020.7500	Label with Plasma indicator	Unit: 7.500 labels
REF: 606.030.7500	Label with EO indicator	Unit: 7.500 labels
REF: 606.040.7500	Label with Form indicator	Unit: 7.500 labels
REF: 606.101.7500	Label blue color & Steam Ind	Unit: 7.500 labels
REF: 606.102.7500	Label red color & Steam Ind	Unit: 7.500 labels
REF: 606.103.7500	Label green color& Steam Ind	Unit: 7.500 labels
REF: 606.104.7500	Label yellow color& Steam Ind	Unit: 7.500 labels
REF: 606.105.7500	Label orange color & Steam Ind	Unit: 7.500 labels
REF: 606.106.7500	Label purple color & Steam Ind	Unit: 7.500 labels

